

MUCOSAL ATOMIZATION DEVICE (MAD)

BASIC EMT / ADVANCED EMT / PARAMEDIC

1. May be used to administer Versed, Ativan, Fentanyl, and Narcan (**based on level of certification**) to adult and pediatric patients.
2. Draw up an extra 0.1 ml of the medication to account for the dead space in the device.
3. Maximum total drug delivery is 2 ml.
4. Briskly compress the syringe plunger to deliver half of the medication into one nostril.
5. Administer the remaining medication into the other nostril.

